

Protocol leukocyte isolation without platelets

(this protocol allows phenotyping cells for tube 1 and 2)

Blood extraction with EDTA

1. Dilute buffy coat **1:5** with sterile PBS (PBS – without Ca and Mg plus **2 mM EDTA**).
2. Add 20 mL of 4% dextran solution to 20 mL of the diluted blood in a 50 mL tube. (2 tubes if 40 mL of blood)
3. Mix and let erythrocytes sediment for 30 min.
4. Collect the supernatant carefully (avoid to take erythrocytes) and centrifuge it at **880 rpm for 15 min (no break)**.

Lysis step

1. Remove the supernatant.
2. Add 7,5 mL of 0,2% NaCl solution to the pellet for 55 seconds (keep mixing the cells with the hypotonic solution using a 10 mL pipette)
3. Add 17,5 mL of 1,2% NaCl solution and then PBS to get a final volume of 50 mL.
4. Centrifuge at **880 rpm for 15 min (break = 0)**
5. Repeat the lysis step once more since there are still plenty of erythrocytes.
6. Centrifuge at **880 rpm for 15 min (break = 0)**
7. Resuspend cells with PBS + 0.09% NaN₃+0.5% BSA+ 1mM EDTA.
8. Count WBC concentration (make a note into protocol Backbone HLDA table)
9. **Centrifuge 8min, 500g, RT, remove the supernatant carefully**
10. Resuspend cells with PBS + 0.09% NaN₃+0.5% BSA (final concentration **40 million/mL**)

Procedure for staining, see the protocol Backbone HLDA table.xls (stainings tab and 96-well tab)

Sample staining (work in 96 well plate)

11. Add Test mAb into each row as per manufacturer recommendation titer (5ul or 10ul, see Backbone HLDA table.xls, 96-well tab)
12. Pipette 40ul of cell ($1,6 \times 10^6$ cells) suspension into V-bottom 96-well plate
13. Where 5ul (or 2.5ul) of test mAb is recommended, add 5ul (or 7.5ul) of PBS+0.09% of NaN₃+0.5% BSA
14. Incubate for 30 min at RT protected from light
15. Prepare reagent mix and add PBS to reach 25ul mix per well (see protocol Backbone HLDA table.xls, Stainings tab)
16. Mix reagent mix with pipette and add 25ul of reagent mix into 50ul of cell suspension, mix with pipette.
17. Incubate for 30 min at RT protected from light
18. Add 100 μ l PBS+0.09% of NaN₃+0.5% BSA
19. Centrifuge 8min, 500g, RT, dump to sink
20. Resuspend the cell pellet in 200 μ l PBS+0.09% of NaN₃+0.5% BSA.
21. Centrifuge 8min, 500g, RT, dump to sink
22. Resuspend the cell pellet in 200 μ l PBS+2mmol EDTA.
23. Reconstitute the Quantibrite PE beads with 250ul PBS+0.09% of NaN₃+0.5% of BSA, vortex and add 200ul to the plate
24. Acquire the cells after staining or (if not immediately acquired) store at darkness for max 1 h until measured in the flow cytometer
25. Acquire tube on HTS using setting listed in protocol Backbone HLDA table.xls, 96-well tab
26. Block lid sensor with magnet, pipette up & down 200ul with multichannel pipette each row just before it is acquired by HTS.

Solutions:

All the solutions are made with sterile bidistilled water using sterile material to avoid the presence of endotoxin. The solutions have to be filtered with 0,2 μm filters.

- 4% dextran solution (in 0.9% NaCl): 0,9g NaCl, 4g dextran in 100 mL bidistilled water.
- 0.2% NaCl solution: 0,2g NaCl in 100 mL bidistilled water.
- 1.2% NaCl solution: 1,2g NaCl in 100 mL bidistilled water.
- PBS without Ca and Mg plus 2 mM EDTA.
- PBS+0.09% of NaN₃+0.5% BSA

Note: Do not decant the tubes. Always use a pipette to discard the supernatants.